



## Safety Data Sheet according to WHS Regulations

Date of issue: 30.06.2025

Revision: 30.06.2025

### 1 Identification

- **Other means of identification**
- **Trade name:** Transcend™
- **Article number:** SDS 441-001.04R01, 4726, 1006864, 1006865, 1006866, 1006867, 1006868, 1006869, 1006870
- **Relevant identified uses of the substance or mixture and uses advised against** Professional Dental Composite
- **Application of the substance / the mixture** Professional Dental Composite
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Ultradent Products, Inc.  
505 W. Ultradent Drive (10200 S)  
South Jordan, UT 84095-3942  
USA  
onlineordersupport@ultradent.com  
(800) 552-5512  
  
Ultradent Australia Pty Ltd.  
Level 22/2 Market Street  
Sydney NSW 2000  
Australia  
Email: info.anz@ultradent.com  
Toll Free: 1-800-290929
- **Further information obtainable from:** Customer Service
- **Emergency telephone number:**  
CHEMTREC (NORTH AMERICA) : +1 (800) 424-9300  
(INTERNATIONAL) : +(703) 527-3887

### 2 Hazard(s) Identification

- **Classification of the substance or mixture**



Skin sensitisation – Category 1 H317 May cause an allergic skin reaction.

- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** GHS07
- **Signal word** Warning
- **Hazard-determining components of labelling:**  
Hydroxypropyl Methacrylate (≥1-<10 %)
- **Hazard statements**  
H317 May cause an allergic skin reaction.
- **Precautionary statements**  
P101 If medical advice is needed, have product container or label at hand.  
P102 Keep out of reach of children.  
P103 Read label before use.  
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.  
P280 Wear protective gloves.  
P362+P364 Take off contaminated clothing and wash it before reuse.

(Contd. on page 2)

# Safety Data Sheet

## according to WHS Regulations

Date of issue: 30.06.2025

Revision: 30.06.2025

Trade name: **Transcend™**

(Contd. of page 1)

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P321 Specific treatment (see on this label).

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

### 3 Composition and Information on Ingredients

- **Chemical characterisation: Mixtures**
- **Description:** Mixture of substances listed below with nonhazardous additions.

- **Dangerous components:**

27813-02-1	Hydroxypropyl Methacrylate ⚠ Eye damage/irritation – Category 2A, H319; Skin sensitisation – Category 1, H317	≥1-<10%
6606-59-3	HDDMA ⚠ Skin corrosion/irritation – Category 2, H315; Specific target organ toxicity (single exposure) – Category 3, H335	≥0-<10%

- **Additional information:**

For the wording of the listed hazard phrases refer to section 16.

Not applicable

None

### 4 First Aid Measures

- **General information:** Immediately remove any clothing soiled by the product.
- **After inhalation:**  
This product is a thick paste, therefore inhalation is extremely unlikely.  
Seek medical treatment in case of complaints.  
Supply fresh air and to be sure call for a doctor.  
In case of unconsciousness place patient stably in side position for transportation.
- **After skin contact:** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** If swallowed in large quantities seek medical attention.
- **Information for doctor:**  
· **Most important symptoms and effects, both acute and delayed** No further relevant information available.  
· **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Fire Fighting Measures

- **Suitable extinguishing agents:**  
Water spray, fog, or mist, foam, dry chemical, carbon dioxide  
Use fire extinguishing methods suitable to surrounding conditions.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Protective equipment:** Wear fully protective suit.

### 6 Accidental Release Measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**  
Dispose contaminated material as waste according to section 13.

(Contd. on page 3)

## Safety Data Sheet according to WHS Regulations

Date of issue: 30.06.2025

Revision: 30.06.2025

**Trade name: Transcend™**

(Contd. of page 2)

- Ensure adequate ventilation.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

### 7 Handling and Storage

- **Handling:**
- **Precautions for safe handling:** Ensure good ventilation/exhaustion at the workplace.
- **Information about fire - and explosion protection:** No special measures required.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**  
Protect from exposure to the light.  
See product labelling.
- **Specific end use(s)** Professional Dental Composite

### 8 Exposure controls and personal protection

- **Appropriate engineering controls** No further data; see section 7.
- **Ingredients with limit values that require monitoring at the workplace:**  
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists valid during the making were used as basis.
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
Immediately remove all soiled and contaminated clothing.  
Wash hands before breaks and at the end of work.
- **Respiratory protection:**  
In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.
- **Protection of hands:**



Protective gloves

- The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.  
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.  
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Material of gloves**  
The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
- **Penetration time of glove material**  
The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Not required.

(Contd. on page 4)

AU

# Safety Data Sheet

## according to WHS Regulations

Date of issue: 30.06.2025

Revision: 30.06.2025

Trade name: Transcend™

(Contd. of page 3)

· **Body protection:** Protective work clothing

### 9 Physical and Chemical Properties

#### · General Information

#### · Appearance:

· <b>Form:</b>	Paste
· <b>Colour:</b>	Shade dependent
· <b>Odour:</b>	Acrylic
· <b>Odour threshold:</b>	Not determined.
· <b>pH-value:</b>	Not applicable (non-aqueous)
· <b>Change in condition</b>	
· <b>Melting point/freezing point:</b>	Undetermined.
· <b>Initial boiling point and boiling range:</b>	Undetermined.
· <b>Flash point:</b>	Not applicable.
· <b>Flammability</b>	Not determined.
· <b>Decomposition temperature:</b>	Not determined.
· <b>Ignition temperature:</b>	Product is not selfigniting.
· <b>Explosive properties:</b>	Product does not present an explosion hazard.
· <b>Explosion limits:</b>	
· <b>Lower:</b>	Not determined.
· <b>Upper:</b>	Not determined.
· <b>Vapour pressure:</b>	Not applicable.
· <b>Density at 20 °C:</b>	2.04 g/cm <sup>3</sup>
· <b>Relative density</b>	Not determined.
· <b>Vapour density</b>	Not applicable.
· <b>Evaporation rate</b>	Not applicable.
· <b>Solubility in / Miscibility with</b>	
· <b>water:</b>	Insoluble.
· <b>Partition coefficient: n-octanol/water:</b>	Not determined.
· <b>Viscosity:</b>	
· <b>Dynamic:</b>	Not applicable.
· <b>Kinematic:</b>	Not applicable.

#### · Other information

· <b>Particle characteristics</b>	Not determined.
· <b>Physical state</b>	Solid

### 10 Stability and Reactivity

- **Reactivity** No further relevant information available.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions:** No dangerous reactions known.
- **Conditions to avoid:**
  - Light
  - Heat
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

AU

(Contd. on page 5)

# Safety Data Sheet

## according to WHS Regulations

Date of issue: 30.06.2025

Revision: 30.06.2025

Trade name: **Transcend™**

(Contd. of page 4)

### 11 Toxicological Information

- **Information on toxicological effects**
- **Acute toxicity** Based on available data, the classification criteria are not met.
- **Primary irritant effect:**
- **Skin corrosion/irritation** Based on available data, the classification criteria are not met.
- **Serious eye damage/irritation** Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitisation** May cause an allergic skin reaction.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

### 12 Ecological Information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behaviour in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:**  
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water  
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**  
Dispose of contents/container in accordance with international, federal, state, and local regulations.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

### 14 Transport information

- |                                  |                      |
|----------------------------------|----------------------|
| · <b>UN-Number</b>               |                      |
| · <b>ADG, IMDG, IATA</b>         | <i>not regulated</i> |
| · <b>UN proper shipping name</b> |                      |
| · <b>ADG, ADN, IMDG, IATA</b>    | <i>not regulated</i> |

(Contd. on page 6)

AU

# Safety Data Sheet

## according to WHS Regulations

Date of issue: 30.06.2025

Revision: 30.06.2025

**Trade name: Transcend™**

(Contd. of page 5)

· <b>Transport hazard class(es)</b>	
· <b>ADG, ADN, IMDG, IATA</b>	
· <b>Class</b>	not regulated
· <b>Packing group</b>	
· <b>ADG, IMDG, IATA</b>	not regulated
· <b>Environmental hazards:</b>	Not applicable.
· <b>Special precautions for user</b>	Not Applicable
· <b>Transport in bulk according to Annex II of Marpol and the IBC Code</b>	Not applicable.
· <b>UN "Model Regulation":</b>	not regulated

### 15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

75-56-9	Methyloxirane
---------	---------------

· **Australian Inventory of Industrial Chemicals**

65997-17-3	Inorganic Glass
41637-38-1	EOBPADMA
	Trade Secret
27813-02-1	Hydroxypropyl Methacrylate
6606-59-3	HDDMA
79-41-4	Methacrylic Acid
13463-67-7	Titanium Dioxide
10287-53-3	Ethyl-4-Dimethylamino Benzoate
	Organophosphine Oxide
75-56-9	Methyloxirane
10373-78-1	Camphorquinone
51274-00-1	Yellow Iron Oxide
	Proprietary
	Chromium Iron Oxide
	Proprietary
	Proprietary
1332-37-2	Red Iron Oxide

· **Standard for the Uniform Scheduling of Medicines and Poisons**

75-56-9	Methyloxirane	S7
	Proprietary	S4, S6

· **Australia: Priority Existing Chemicals**

None of the ingredients is listed.

· **Directive 2012/18/EU**

· **Named dangerous substances - ANNEX I** None of the ingredients is listed.

(Contd. on page 7)

# Safety Data Sheet

## according to WHS Regulations

Date of issue: 30.06.2025

Revision: 30.06.2025

**Trade name: Transcend™**

(Contd. of page 6)

· **Chemical safety assessment:**

Device is biocompatible when used as directed by dental professionals per ISO 10993-1

### 16 Other information

*This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.*

· **Relevant phrases from Section 3**

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

H335 May cause respiratory irritation.

· **Department issuing SDS:** Environmental, Health, and Safety

· **Contact:** Customer Service

· **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

Skin corrosion/irritation – Category 2: Skin corrosion/irritation – Category 2

Eye damage/irritation – Category 2A: Serious eye damage/eye irritation – Category 2A

Skin sensitisation – Category 1: Skin sensitisation – Category 1

Specific target organ toxicity (single exposure) – Category 3: Specific target organ toxicity (single exposure) – Category 3

· **\* Data compared to the previous version altered.**

AU